

For Research Use Only

# Anti-Human STRO-1 (STRO-1)

Catalog Number: 65184-1-Ig



## Basic Information

Catalog Number:

65184-1-Ig

Size:

100ug, 0.5 mg/ml

Source:

Mouse

Isotype:

IgM, lambda

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Purified by precipitation and chromatography

CloneNo.:

STRO-1

## Applications

Tested Applications:

FC

Species Specificity:

Human

## Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> lysed human blood cells were surface stained with 0.2 ug Anti-Human STRO-1 (65184-1-Ig, Clone:STRO-1) and AlexaFluor-Goat anti-Mouse IgG at dilution 1:1000. Cells were not fixed. CD45 dim, CD34 positive cells were gated.



1X10<sup>6</sup> lysed human blood cells were surface stained with 0.2 ug Anti-Human STRO-1 (65184-1-Ig, Clone:STRO-1) and AlexaFluor-Goat anti-Mouse IgG at dilution 1:1000. Cells were not fixed. CD45 dim, CD34 positive cells were gated.