

For Research Use Only

Anti-Mouse Integrin beta-7 (FIB504)

Catalog Number: 65213-1-Ig



Basic Information

Catalog Number: 65213-1-Ig	GenBank Accession Number: BC011184	Purification Method: Protein A purification
Size: 100ug, 0.5 mg/ml	GeneID (NCBI): 16421	CloneNo.: FIB504
Source: Rat	UNIPROT ID: P26011	
Isotype: IgG2a, kappa	Full Name: integrin beta 7	

Applications

Tested Applications:
FC

Species Specificity:
Human, Mouse

Storage


Storage:
Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.1% sodium azide, pH 7.3.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ C57BL/6 mouse splenocytes were surface stained with 0.5 ug Anti-Mouse Integrin beta-7 (65213-1-Ig, Clone: FIB504) or 0.5 ug Rat IgG2a Isotype Control (2A3) (65209-1-Ig, Clone: 2A3) and FITC anti-Rat IgG2a Antibody at dilution 1:100. Cells were co-stained with CoraLite® Plus 647 Anti-Mouse CD3 (17A2) (CL647-65077, Clone: 17A2). Cells were not fixed.