

For Research Use Only

Anti-Human CD96 (NK92.39)

Catalog Number: 65250-1-Ig



Basic Information

Catalog Number:

65250-1-Ig

Size:

100ug, 0.5 mg/ml

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

BC020749

GeneID (NCBI):

10225

UNIPROT ID:

P40200

Full Name:

CD96 molecule

Calculated MW:

402 aa, 46 kDa

Purification Method:

Affinity purification

CloneNo.:

NK92.39

Applications

Tested Applications:

FC

Species Specificity:

Human

Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human PBMCs were surface co-stained with CL750 Anti-Human CD19 and 0.2 ug Anti-Human CD96 (65250-1-Ig, Clone: NK92.39) and FITC-Donkey anti-Mouse IgG at dilution 1:1000, or Mouse IgG1 Isotype Control. Cells were not fixed. Lymphocytes were gated.



1X10⁶ human PBMCs were surface stained with 0.2 ug Anti-Human CD96 (65250-1-Ig, Clone: NK92.39) and FITC-Donkey anti-Mouse IgG at dilution 1:1000, or Mouse IgG1 Isotype Control. Cells were not fixed. Lymphocytes were gated.