

For Research Use Only

Anti-Mouse GL7 (GL-7)

Catalog Number: 65261-1-Ig



Basic Information

Catalog Number:

65261-1-Ig

Size:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgM, kappa

GenBank Accession Number:

GeneID (NCBI):

104230

Full Name:

lymphocyte antigen 77

Purification Method:

Affinity purification

CloneNo.:

GL-7

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ unstimulated (blue) or plate-bound anti-CD3 and anti-CD28 stimulated (red) mouse splenocytes were surface stained with 0.5 ug Anti-Mouse GL7 (65261-1-Ig, Clone:GL-7) and FITC anti-rat IgM Antibody at dilution 1:100. Cells were not fixed.