

For Research Use Only

Anti-Rat CD43 (W3/13) Mouse IgG2a Recombinant Antibody, PBS Only

Catalog Number: 65607-1-PBS



Basic Information

Catalog Number:

65607-1-PBS

Size:

1mg, 2 mg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

GeneID (NCBI):

24796

Full Name:

sialophorin

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

W3/13

Applications

Tested Applications:

FC

Species Specificity:

rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ rat splenocytes were surface stained with 0.25 ug Anti-Rat CD43 Mouse IgG2a Recombinant Antibody (65607-1-MR, Clone: W3/13) or Mouse IgG2a Isotype Control (65208-1-Ig, Clone: C1.18.4) and FITC anti-Mouse IgG2a Antibody. Cells were not fixed. This data was developed using the same antibody clone with 65607-1-PBS in a different storage buffer formulation.



1x10⁶ rat splenocytes were surface stained with 0.25 ug Anti-Rat CD43 Mouse IgG2a Recombinant Antibody (65607-1-MR, Clone: W3/13) (red) or Mouse IgG2a Isotype Control (65208-1-Ig, Clone: C1.18.4) (blue) and FITC anti-Mouse IgG2a Antibody. Cells were not fixed. This data was developed using the same antibody clone with 65607-1-PBS in a different storage buffer formulation.