

For Research Use Only

Anti-Human CD1c (L161) Mouse IgG2a Recombinant Antibody, PBS Only

Catalog Number: 65626-1-PBS



Basic Information

Catalog Number:

65626-1-PBS

Size:

1mg, 2mg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

BC126465

GeneID (NCBI):

911

Full Name:

CD1c molecule

Purification Method:

Protein A purification

CloneNo.:

L161

Applications

Tested Applications:

FC

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human peripheral blood lymphocytes were surface stained with 0.25 ug Anti-Human CD1c (L161) Mouse IgG2a Recombinant Antibody (65626-1-MR, Clone: L161) or Mouse IgG2a Isotype Control (C1.18.4) (65208-1-Ig, Clone: C1.18.4) and PE-conjugated Goat Anti-Mouse IgG. Cells were co-stained with CoraLite® Plus 647 Anti-Human CD19. Cells were not fixed. This data was developed using the same antibody clone with 65626-1-PBS in a

