

For Research Use Only

NDUFA4L2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66050-5-PBS



Basic Information

Catalog Number: 66050-5-PBS	GenBank Accession Number: BC011910	Purification Method: Protein A Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 56901	CloneNo.: 3C1C5
Source: Mouse	UNIPROT ID: Q9NRX3	
Isotype: IgG2b	Full Name: NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 4-like 2	
Immunogen Catalog Number: AG9233	Calculated MW: 87 aa, 10 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

66050-5-PBS targets NDUFA4L2 as part of a matched antibody pair:

MP50529-3: 66050-4-PBS capture and 66050-5-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50529-3, NDUF4L2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66050-4-PBS. Detection antibody: 66050-5-PBS. Standard: Ag9233. Range: 0.195-100 ng/mL