

For Research Use Only

MBIP Monoclonal antibody, PBS Only



Catalog Number: 66102-1-PBS

Featured Product

Basic Information

Catalog Number: 66102-1-PBS	GenBank Accession Number: BC005197	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 51562	CloneNo.: 2E6H11
Source: Mouse	UNIPROT ID: Q9NS73	
Isotype: IgG2b	Full Name: MAP3K12 binding inhibitory protein 1	
Immunogen Catalog Number: AG18090	Calculated MW: 39 kDa	
	Observed MW: 39 kDa	

Applications

Tested Applications:
WB, IP, IHC, Indirect ELISA

Species Specificity:
human, mouse

Background Information

MBIP (MAP3K12 binding inhibitory protein 1), also known as MAPK upstream kinase (MUK)-binding inhibitory protein, is a 344 amino acid nuclear and cytoplasmic protein that assists in the activation of the JNK pathway by inhibiting ZPK activity. Ubiquitously expressed, MBIP is found at high levels in lung and heart and is one component of the ADA2A-containing complex (ATAC) in conjunction with ADA3, GCN5, CDC101, TADA2L, WDR5, YEATS2, DR1 and CRP2BP.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 66102-1-Ig (MBIP antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66102-1-PBS in a different storage buffer formulation.




WB result of MBIP antibody (66102-1-Ig; 1:2000; room temperature for 1.5 hours) with sh-Control and sh-MBIP transfected HEK-293 cells. This data was developed using the same antibody clone with 66102-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 66102-1-Ig (MBIP antibody) at dilution of 1:800 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 66102-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 66102-1-Ig (MBIP antibody) at dilution of 1:800 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 66102-1-PBS in a different storage buffer formulation.



IP result of anti-MBIP (IP:66102-1-Ig, 3ug; Detection:66102-1-Ig 1:300) with HEK-293 cells lysate 1800ug. This data was developed using the same antibody clone with 66102-1-PBS in a different storage buffer formulation.