

For Research Use Only

# CD29 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66315-2-PBS

## Basic Information

**Catalog Number:**

66315-2-PBS

**Size:**

100ug, Concentration: 1 mg/ml by  
Nanodrop;

**Source:**

Mouse

**Isotype:**

IgG2a

**Immunogen Catalog Number:**

AG24476

**GenBank Accession Number:**

BC020057

**GeneID (NCBI):**

3688

**ENSEMBL Gene ID:**

ENSG00000150093

**UNIPROT ID:**

P05556

**Full Name:**

integrin, beta 1 (fibronectin receptor,  
beta polypeptide, antigen CD29  
includes MDF2, MSK12)

**Calculated MW:**

88 kDa

**Purification Method:**

Protein A purification

**CloneNo.:**

2A1F6

## Applications

**Tested Applications:**

Cytometric bead array, Indirect ELISA

**Species Specificity:**

human

## Storage

**Storage:**

Store at -80°C.

**Storage Buffer:**

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50343-1, CD29 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66315-2-PBS. Detection antibody: 66315-3-PBS. Standard: Ag24476. Range: 3.125-200 ng/mL.