

For Research Use Only

TrkC Monoclonal antibody, PBS Only

Catalog Number: 66380-1-PBS



Basic Information

Catalog Number:

66380-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG2635

GenBank Accession Number:

BC013693

GeneID (NCBI):

4916

UNIPROT ID:

Q16288

Full Name:

neurotrophic tyrosine kinase, receptor, type 3

Calculated MW:

839 aa, 94 kDa

Observed MW:

110 kDa, 160 kDa

Purification Method:

Protein G purification

CloneNo.:

2A7F11

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, rat, mouse, pig

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66380-1-Ig (TrkC antibody at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66380-1-PBS in a different storage buffer formulation.