

For Research Use Only

# OMA1 Monoclonal antibody, PBS Only (Capture)

Catalog Number:66572-1-PBS



## Basic Information

<b>Catalog Number:</b> 66572-1-PBS	<b>GenBank Accession Number:</b> BC012915	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 100ug , Concentration: 1mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 115209	<b>CloneNo.:</b> 2E11D7
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q96E52	
<b>Isotype:</b> IgG1	<b>Full Name:</b> OMA1 homolog, zinc metallopeptidase (S. cerevisiae)	
<b>Immunogen Catalog Number:</b> AG11201	<b>Calculated MW:</b> 524 aa, 60 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

66572-1-PBS targets OMA1 as part of a matched antibody pair:

MP50571-1: 66572-1-PBS capture and 66572-2-PBS detection (validated in Cytometric bead array)

MP50571-3: 66572-1-PBS capture and 67449-1-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50571-1, OMA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66572-1-PBS. Detection antibody: 66572-2-PBS. Standard:Ag11201. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP50571-3, OMA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66572-1-PBS. Detection antibody: 67449-1-PBS. Standard:Ag11201. Range: 3.125-100 ng/mL