

For Research Use Only

CHCHD6 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66597-3-PBS



Basic Information

Catalog Number: 66597-3-PBS	GenBank Accession Number: BC006123	Purification Method: Protein G Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 84303	CloneNo.: 2G10C12
Source: Mouse	UNIPROT ID: Q9BRQ6	
Isotype: IgG1	Full Name: coiled-coil-helix-coiled-coil-helix domain containing 6	
Immunogen Catalog Number: AG14700	Calculated MW: 235 aa, 26 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

66597-3-PBS targets CHCHD6 as part of a matched antibody pair:

MP50477-1: 66597-2-PBS capture and 66597-3-PBS detection (validated in Cytometric bead array)

MP50477-4: 66597-1-PBS capture and 66597-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50477-1, CHCHD6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66597-2-PBS. Detection antibody: 66597-3-PBS. Standard: Ag14700. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50477-4, CHCHD6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66597-1-PBS. Detection antibody: 66597-3-PBS. Standard: Ag14700. Range: 0.391-100 ng/mL.