

For Research Use Only

GRIK2 Monoclonal antibody, PBS Only



Catalog Number: 66631-2-PBS

Basic Information

Catalog Number:

66631-2-PBS

Size:

100ug, Concentration: 1mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25308

GenBank Accession Number:

BC037954

GeneID (NCBI):

2898

UNIPROT ID:

Q13002

Full Name:

glutamate receptor, ionotropic, kainate 2

Calculated MW:

908 aa, 103 kDa

Observed MW:

103-115 kDa

Purification Method:

Protein G purification

CloneNo.:

1E5G6

Applications

Tested Applications:

Indirect ELISA, WB

Species Specificity:

pig, rat, mouse, human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66631-2-Ig (GRIK2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66631-2-PBS in a different storage buffer formulation.