

For Research Use Only

CD302 Monoclonal antibody, PBS Only



Catalog Number: 66640-1-PBS

Basic Information

Catalog Number: 66640-1-PBS	GenBank Accession Number: BC020646	Purification Method: Protein G purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 9936	CloneNo.: 2B1E8
Source: Mouse	UNIPROT ID: Q8IX05	
Isotype: IgG1	Full Name: CD302 molecule	
Immunogen Catalog Number: AG9873	Calculated MW: 170aa,19 kDa; 232aa,26 kDa	
	Observed MW: 30-32 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, Mouse

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



THP-1 cells were subjected to SDS PAGE followed by western blot with 66640-1-Ig (CD302 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66640-1-PBS in a different storage buffer formulation.