

For Research Use Only

CD133 Monoclonal antibody, PBS Only (Capture)

Catalog Number:66666-2-PBS



Basic Information

Catalog Number: 66666-2-PBS	GenBank Accession Number: BC012089	Purification Method: Protein G purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 8842	CloneNo.: 3A6B8
Source: Mouse	UNIPROT ID: O43490	
Isotype: IgG1	Full Name: prominin 1	
Immunogen Catalog Number: AG13327	Calculated MW: 97 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

66666-2-PBS targets CD133 as part of a matched antibody pair:

MP50457-1: 66666-2-PBS capture and 66666-3-PBS detection (validated in Cytometric bead array)

MP50457-2: 66666-2-PBS capture and 66666-4-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50457-2, CD133 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66666-2-PBS. Detection antibody: 66666-4-PBS. Standard:Ag13327. Range: 3.125-100 ng/mL

Cytometric bead array standard curve of MP50457-1, CD133 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66666-2-PBS. Detection antibody: 66666-3-PBS. Standard:Ag13327. Range: 1.563-100 ng/mL