

For Research Use Only

ABCD3 Monoclonal antibody, PBS Only



Catalog Number: 66697-1-PBS

Basic Information

Catalog Number: 66697-1-PBS	GenBank Accession Number: BC068509	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 5825	CloneNo.: 1G7G5
Source: Mouse	UNIPROT ID: P28288	
Isotype: IgG3	Full Name: ATP-binding cassette, sub-family D (ALD), member 3	
Immunogen Catalog Number: AG24760	Observed MW: 70 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66697-1-Ig (ABCD3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66697-1-PBS in a different storage buffer formulation.