

For Research Use Only

CD21 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66701-6-PBS



Basic Information

Catalog Number: 66701-6-PBS	GenBank Accession Number: BC136394	Purification Method: Protein A Magarose purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 1380	CloneNo.: 1H7A4
Source: Mouse	ENSEMBL Gene ID: ENSG00000117322	
Isotype: IgG2a	UNIPROT ID: P20023	
Immunogen Catalog Number: AG18058	Full Name: complement component (3d/Epstein Barr virus) receptor 2	
	Calculated MW: 1092 aa, 119 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

66701-6-PBS targets CD21 as part of a matched antibody pair:

MP51164-3: 66701-2-PBS capture and 66701-6-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51164-3, CD21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66701-2-PBS. Detection antibody: 66701-6-PBS. Standard: Ag18058. Range: 0.781-100 ng/mL.