

For Research Use Only

PRDX3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66810-3-PBS



Basic Information

Catalog Number: 66810-3-PBS	GenBank Accession Number: BC007062	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 10935	CloneNo.: 1F1A9
Source: Mouse	UNIPROT ID: P30048	
Isotype: IgG2b	Full Name: peroxiredoxin 3	
Immunogen Catalog Number: AG1062	Calculated MW: 27 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

66810-3-PBS targets PRDX3 as part of a matched antibody pair.

MP50218-1: 66810-2-PBS capture and 66810-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50218-1, PRDX3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66810-2-PBS. Detection antibody: 66810-3-PBS. Standard: Ag1062. Range: 0.781-100 ng/mL.