

For Research Use Only

IBA1 Monoclonal antibody

Catalog Number: 66827-2-Ig **1 Publications**



Basic Information

Catalog Number: 66827-2-Ig	GenBank Accession Number: BC009474	Purification Method: Protein G purification
Size: 150ul , Concentration: 1000 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 199	CloneNo.: 2B2A12
Source: Mouse	ENSEMBL Gene ID: ENSG00000204472	Recommended Dilutions: IHC 1:2000-1:8000
Isotype: IgG1	Full Name: allograft inflammatory factor 1	
Immunogen Catalog Number: AG28236	Calculated MW: 17 kDa	

Applications

Tested Applications: IHC, ELISA	Positive Controls: IHC : human brain tissue, human liver cancer tissue, human tonsillitis tissue
Cited Applications: IF	
Species Specificity: Human	
Cited Species: human	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Notable Publications

Author	Pubmed ID	Journal	Application
Sara Serafini	38397814	Antioxidants (Basel)	IF

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data