

For Research Use Only

Decorin Monoclonal antibody, PBS Only



Catalog Number: 66847-1-PBS

Basic Information

Catalog Number: 66847-1-PBS	GenBank Accession Number: BC005322	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 1634	CloneNo.: 3A4E4
Source: Mouse	UNIPROT ID: P07585	
Isotype: IgG1	Full Name: decorin	
Immunogen Catalog Number: AG6275	Calculated MW: 359 aa, 40 kDa	
	Observed MW: 40 kDa, 100 kDa	

Applications

Tested Applications:
WB, IF, ELISA

Species Specificity:
Human, Pig

Background Information

Decorin is a member of the small leucine-rich proteoglycan family of proteins. Decorin interacts with type I collagen fibrils, thereby influencing the kinetics of fibril formation and the distance between adjacent collagen fibrils. Binding of this protein to multiple cell surface receptors mediates its role in tumor suppression, including a stimulatory effect on autophagy and inflammation and an inhibitory effect on angiogenesis and tumorigenesis.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



pig cartilage tissue were subjected to SDS PAGE followed by western blot with 66847-1-Ig (Decorin antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66847-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed BxPC-3 cells using Decorin antibody (66847-1-Ig, Clone: 3A4E4) at dilution of 1:800 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 66847-1-PBS in a different storage buffer formulation.