

For Research Use Only

# Integrin alpha 6 Monoclonal antibody, PBS Only



Catalog Number: 66906-1-PBS

## Basic Information

<b>Catalog Number:</b> 66906-1-PBS	<b>GenBank Accession Number:</b> BC136456	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug, Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 3655	<b>CloneNo.:</b> 2B12A1
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> P23229	
<b>Isotype:</b> IgG2b	<b>Full Name:</b> integrin, alpha 6	
<b>Immunogen Catalog Number:</b> AG25866	<b>Calculated MW:</b> 1073 aa, 119 kDa	
	<b>Observed MW:</b> 115-120 kDa	

## Applications

**Tested Applications:**  
WB, Indirect ELISA

**Species Specificity:**  
Human

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 66906-1-Ig (Integrin alpha-6 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66906-1-PBS in a different storage buffer formulation.