

For Research Use Only

Integrin alpha-6 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number:66906-2-PBS



Basic Information

Catalog Number: 66906-2-PBS	GenBank Accession Number: BC136456	Purification Method: Protein A Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 3655	CloneNo.: 2G2D4
Source: Mouse	UNIPROT ID: P23229	
Isotype: IgG2a	Full Name: integrin, alpha 6	
Immunogen Catalog Number: AG25866	Calculated MW: 1073 aa, 119 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

66906-2-PBS targets Integrin alpha-6 as part of a matched antibody pair:

MP50396-1: 66906-2-PBS capture and 66906-3-PBS detection (validated in Cytometric bead array)

MP50396-2: 66906-4-PBS capture and 66906-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50396-1, Integrin alpha-6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66906-2-PBS. Detection antibody: 66906-3-PBS. Standard: Ag25866. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50396-2, Integrin alpha-6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66906-4-PBS. Detection antibody: 66906-2-PBS. Standard: Ag25866. Range: 3.125-100 ng/mL.