

For Research Use Only

NDC80 Monoclonal antibody, PBS Only



Catalog Number: 66960-1-PBS

Basic Information

Catalog Number: 66960-1-PBS	GenBank Accession Number: BC035617	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 10403	CloneNo.: 2B7A8
Source: Mouse	UNIPROT ID: O14777	
Isotype: IgG1	Full Name: NDC80 homolog, kinetochore complex component (S. cerevisiae)	
Immunogen Catalog Number: AG8065	Calculated MW: 74 kDa	
	Observed MW: 71-74 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, rat, mouse

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66960-1-Ig (NDC80 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66960-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 66960-1-Ig (NDC80 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66960-1-PBS in a different storage buffer formulation.