

For Research Use Only

EPAC2 Monoclonal antibody, PBS Only



Catalog Number: 67044-1-PBS

Basic Information

Catalog Number: 67044-1-PBS	GenBank Accession Number: BC024004	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 11069	CloneNo.: 1H4A10
Source: Mouse	UNIPROT ID: Q8WZA2	
Isotype: IgG1	Full Name: Rap guanine nucleotide exchange factor (GEF) 4	
Immunogen Catalog Number: AG5186	Calculated MW: 116 kDa	
	Observed MW: 110 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, mouse, rat

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67044-1-Ig (EPAC2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67044-1-PBS in a different storage buffer formulation.