

For Research Use Only

Tissue Factor Monoclonal antibody, PBS Only



Catalog Number: 67056-1-PBS

Basic Information

Catalog Number: 67056-1-PBS	GenBank Accession Number: BC011029	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 2152	CloneNo.: 2E4A1
Source: Mouse	UNIPROT ID: P13726	
Isotype: IgG1	Full Name: coagulation factor III (thromboplastin, tissue factor)	
Immunogen Catalog Number: AG27563	Calculated MW: 295 aa, 33 kDa	
	Observed MW: 45-50 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 67056-1-Ig (Tissue Factor antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67056-1-PBS in a different storage buffer formulation.