

For Research Use Only

PPP1CA Monoclonal antibody, PBS Only

Catalog Number: 67070-1-PBS



Basic Information

Catalog Number:

67070-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28451

GenBank Accession Number:

BC001888

GeneID (NCBI):

5499

UNIPROT ID:

P62136

Full Name:

protein phosphatase 1, catalytic subunit, alpha isoform

Calculated MW:

330 aa, 38 kDa

Observed MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

1C10A2

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67070-1-Ig (PPP1CA antibody) at dilution of 1:2600 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Alpha Tubulin Monoclonal antibody (HRP-66031) as loading control. This data was developed using the same antibody clone with 67070-1-PBS in a different storage buffer formulation.



zebrafish tissue were subjected to SDS PAGE followed by western blot with 67070-1-Ig (PPP1CA antibody) at dilution of 1:2600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67070-1-PBS in a different storage buffer formulation.