

For Research Use Only

ILT-3/LILRB4 Monoclonal antibody, PBS Only

Catalog Number: 67119-1-PBS



Basic Information

Catalog Number: 67119-1-PBS	GenBank Accession Number: BC026309	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 11006	CloneNo.: 1C8E1
Source: Mouse	UNIPROT ID: Q8NHJ6	
Isotype: IgG2b	Full Name: Leukocyte immunoglobulin-like receptor, subfamily B (with TM and ITIM domains), member 4	
Immunogen Catalog Number: AG3219	Calculated MW: 448 aa, 49 kDa	
	Observed MW: 49 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MCF-7 cells were subjected to SDS PAGE followed by western blot with 67119-1-Ig (LILRB4 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67119-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed A549 cells using LILRB4 antibody (67119-1-Ig, Clone: 1C8E1) at dilution of 1:400 and Multi-rAb CoraLite® Plus 488-Goat Anti-Mouse Recombinant Secondary Antibody (H+L) (RGAM002). This data was developed using the same antibody clone with 67119-1-PBS in a different storage buffer formulation.