

For Research Use Only

Mouse mFC tag Monoclonal antibody, PBS Only



Catalog Number: 67122-1-PBS

Basic Information

Catalog Number:

67122-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27844

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Calculated MW:

26 kDa

Observed MW:

25-30 kDa

Purification Method:

Protein G purification

CloneNo.:

2F5A2

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

mouse

Background Information

This antibody detects the mouse Fc tag. It's mouse Fc gamma 2a.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



67122-1-Ig was lysed and subjected to SDS-PAGE followed by western blot with 67122-1-Ig and isotype control antibody 66360-1-Ig at dilution of 1:5000 incubated at room temperature for 1.5 hours. HRP-goat anti-mouse kappa chain secondary antibody was used for detection. . This data was developed using the same antibody clone with 67122-1-PBS in a different storage buffer formulation.