

For Research Use Only

IRAK4 Monoclonal antibody, PBS Only



Catalog Number: 67180-1-PBS

Basic Information

Catalog Number: 67180-1-PBS	GenBank Accession Number: BC013316	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 51135	CloneNo.: 1A4C4
Source: Mouse	UNIPROT ID: Q9NWZ3	
Isotype: IgG2a	Full Name: interleukin-1 receptor-associated kinase 4	
Immunogen Catalog Number: AG29174	Calculated MW: 52 kDa	
	Observed MW: 52 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67180-1-Ig (IRAK4 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67180-1-PBS in a different storage buffer formulation.