

For Research Use Only

# DTX2 Monoclonal antibody, PBS Only

Catalog Number: 67209-1-PBS



## Basic Information

Catalog Number:

67209-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG12022

GenBank Accession Number:

BC064535

GeneID (NCBI):

113878

UNIPROT ID:

Q86UW9

Full Name:

deltex homolog 2 (Drosophila)

Calculated MW:

253aa, 28 kDa; 575aa, 64 kDa

Observed MW:

68 kDa

Purification Method:

Protein A purification

CloneNo.:

2D4F8

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67209-1-Ig (DTX2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67209-1-PBS in a different storage buffer formulation.