

For Research Use Only

TMEM230 Monoclonal antibody, PBS Only



Catalog Number: 67247-1-PBS

Basic Information

Catalog Number:

67247-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG15080

GenBank Accession Number:

BC011990

GeneID (NCBI):

29058

UNIPROT ID:

Q96A57

Full Name:

chromosome 20 open reading frame 30

Calculated MW:

183 aa, 20 kDa

Observed MW:

16 kDa

Purification Method:

Protein A purification

CloneNo.:

4D3A7

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, Mouse, Rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67247-1-Ig (TMEM230 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67247-1-PBS in a different storage buffer formulation.