

For Research Use Only

ARF5 Monoclonal antibody, PBS Only

Catalog Number: 67279-1-PBS



Basic Information

Catalog Number: 67279-1-PBS	GenBank Accession Number: BC003043	Purification Method: Protein G purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 381	CloneNo.: 2B3B1
Source: Mouse	UNIPROT ID: P84085	
Isotype: IgG1	Full Name: ADP-ribosylation factor 5	
Immunogen Catalog Number: AG7849	Calculated MW: 21 kDa	
	Observed MW: 18-20 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, mouse, rat, pig

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67279-1-Ig (ARF5 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67279-1-PBS in a different storage buffer formulation.