

For Research Use Only

DLX3 Monoclonal antibody, PBS Only

Catalog Number:67282-1-PBS



Basic Information

Catalog Number: 67282-1-PBS	GenBank Accession Number: BC028970	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 1747	CloneNo.: 1H8E1
Source: Mouse	UNIPROT ID: O60479	
Isotype: IgG2a	Full Name: distal-less homeobox 3	
Immunogen Catalog Number: AG29313	Calculated MW: 287 aa, 32 kDa	
	Observed MW: 50 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat, pig

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67282-1-Ig (DLX3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67282-1-PBS in a different storage buffer formulation.