

For Research Use Only

WDR19 Monoclonal antibody, PBS Only



Catalog Number: 67330-1-PBS

Basic Information

Catalog Number: 67330-1-PBS	GenBank Accession Number: BC032578	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 57728	CloneNo.: 2D1H4
Source: Mouse	UNIPROT ID: Q8NEZ3	
Isotype: IgG2b	Full Name: WD repeat domain 19	
Immunogen Catalog Number: AG4599	Calculated MW: 1342 aa, 152 kDa	
	Observed MW: 145 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67330-1-Ig (WDR19 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67330-1-PBS in a different storage buffer formulation.