

For Research Use Only

N-PAC Monoclonal antibody, PBS Only



Catalog Number: 67340-1-PBS

Basic Information

Catalog Number:

67340-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG7133

GenBank Accession Number:

BC064940

GeneID (NCBI):

84656

UNIPROT ID:

Q49A26

Full Name:

cytokine-like nuclear factor n-pac

Calculated MW:

61 kDa

Observed MW:

60 kDa

Purification Method:

Protein A purification

CloneNo.:

1A5E1

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, Mouse, Rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67340-1-Ig (N-PAC antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67340-1-PBS in a different storage buffer formulation.