

For Research Use Only

MUS81 Monoclonal antibody, PBS Only



Catalog Number: 67351-1-PBS

Featured Product

Basic Information

Catalog Number:

67351-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG29089

GenBank Accession Number:

BC009999

GeneID (NCBI):

80198

UNIPROT ID:

Q96NY9

Full Name:

MUS81 endonuclease homolog (S. cerevisiae)

Calculated MW:

61 kDa

Observed MW:

52-72 kDa

Purification Method:

Protein A purification

CloneNo.:

2D9B7

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, Mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



WB result of MUS81 antibody (67351-1-Ig; 1:10000; incubated at room temperature for 1.5 hours) with sh-Control and sh-MUS81 transfected MCF-7 cells. This data was developed using the same antibody clone with 67351-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 67351-1-Ig (MUS81 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67351-1-PBS in a different storage buffer formulation.