

For Research Use Only

CEP164 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67370-2-PBS



Basic Information

Catalog Number: 67370-2-PBS	GenBank Accession Number: BC000602	Purification Method: Protein G purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 22897	CloneNo.: 3B7F2
Source: Mouse	UNIPROT ID: Q9UPV0	
Isotype: IgG1	Full Name: centrosomal protein 164kDa	
Immunogen Catalog Number: AG17699	Calculated MW: 1460 aa, 164 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

67370-2-PBS targets CEP164 as part of a matched antibody pair:

MP50474-1: 67370-2-PBS capture and 67370-3-PBS detection (validated in Cytometric bead array)

MP50474-2: 67370-2-PBS capture and 67370-4-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50474-1, CEP164 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67370-2-PBS. Detection antibody: 67370-3-PBS. Standard: Ag17699. Range: 0.195-50 ng/mL



Cytometric bead array standard curve of MP50474-2, CEP164 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67370-2-PBS. Detection antibody: 67370-4-PBS. Standard: null. Range: null.