

For Research Use Only

SPATA19 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67486-1-PBS



Basic Information

Catalog Number: 67486-1-PBS	GenBank Accession Number: BC058039	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 219938	CloneNo.: 1G12A11
Source: Mouse	UNIPROT ID: Q7Z5L4	
Isotype: IgG1	Full Name: spermatogenesis associated 19	
Immunogen Catalog Number: AG9808	Calculated MW: 167 aa, 19 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

67486-1-PBS targets SPATA19 as part of a matched antibody pair:

MP50970-1: 67486-1-PBS capture and 67486-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50970-1, SPATA19 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67486-1-PBS. Detection antibody: 67486-2-PBS. Standard:Ag9808. Range: 6.25-200 ng/mL