

For Research Use Only

MRE11A Monoclonal antibody, PBS Only



Catalog Number: 67570-1-PBS

Basic Information

Catalog Number: 67570-1-PBS	GenBank Accession Number: BC063458	Purification Method: Protein G purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 4361	CloneNo.: 1B12B2
Source: Mouse	UNIPROT ID: P49959	
Isotype: IgG1	Full Name: MRE11 meiotic recombination 11 homolog A (S. cerevisiae)	
Immunogen Catalog Number: AG28847	Observed MW: 81 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, Mouse, Rat

Background Information

Mre11 (meiotic recombination 11), also named as Mre11A, is a component of the MRN (MRE11/RAD50/NBS1) complex which is a versatile complex, playing a key role in the sensing, processing, and repair of DSBs. MRE11 possesses both endonuclease and 3'-5' exonuclease activity. It contains two DNA-binding domains (DBDs), enabling it to bind both single-stranded DNA (ssDNA) and double-stranded DNA (dsDNA). After DSBs, MRE11 is loaded onto DSBs sites and cleaves DNA by cooperating with CtIP (CtBP [C terminal binding protein] interacting protein) to initiate end resection. It's reported that MRE11 is over-expressed in breast cancers. (PMID:38128537; 22914783)

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67570-1-Ig (MRE11A antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67570-1-PBS in a different storage buffer formulation.