

For Research Use Only

RPS10 Monoclonal antibody, PBS Only



Catalog Number: 67609-1-PBS

Basic Information

Catalog Number: 67609-1-PBS	GenBank Accession Number: BC001032	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 6204	CloneNo.: 1G1E1
Source: Mouse	UNIPROT ID: P46783	
Isotype: IgG2b	Full Name: ribosomal protein S10	
Immunogen Catalog Number: AG7018	Calculated MW: 19 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, mouse, rat


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67609-1-Ig (RPS10 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67609-1-PBS in a different storage buffer formulation.