

For Research Use Only

Virus SARS-CoV-2 Nucleocapsid Phosphoprotein Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67666-1-PBS



Basic Information

Catalog Number: 67666-1-PBS	GenBank Accession Number: NC_045512	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 43740575	CloneNo.: 1B3C3
Source: Mouse	Full Name: COVID-19 N Protein	
Isotype: IgG1		
Immunogen Catalog Number: AG30676		

Applications

Tested Applications:
WB, ELISA, Sandwich ELISA, Indirect ELISA

Species Specificity:
virus

Product Information

67666-1-PBS targets SARS-CoV-2 Nucleocapsid Phosphoprotein as part of a matched antibody pair:

MP50061-1: 67666-1-PBS capture and 67666-2-PBS detection (validated in Sandwich ELISA)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Background Information

The nucleocapsid (N) protein has multiple functions including formation of nucleocapsids, signal transduction virus budding, RNA replication, and mRNA transcription. N protein is an important antigen for coronavirus, and it is normally highly conserved, with a molecular weight of about 50 kDa. it can be used as a marker in diagnostic assays due to its high immunogenicity (PMID: 32416961, PMID: 32235387). 67666-1-Ig can be used as capture antibody. 67666-2-Ig can be used as detection antibody.

Storage




Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Indirect ELISA was carried out by coating eukaryotic expressed N protein at 70 ng/well followed by blocking and adding serial diluted primary antibody 67666-1-Ig and 67666-2-Ig respectively. Signal was developed with TMB and stopped by H2SO4. Signal strength was measured by absorbance at 450 nm. This data was developed using the same antibody clone with 67666-1-PBS in a different storage buffer formulation.

Sandwich ELISA was carried out by coating 67666-1-Ig at 80 ng/well followed by blocking and adding different concentration of eukaryotic expressed N protein (0.5-1000 ng/mL). HRP-conjugated clone 67666-2-Ig was used at 1 µg/mL for detection. Signal was developed with TMB and stopped by H2SO4. Signal strength was measured by absorbance at 450 nm. This data was developed using the same antibody clone with 67666-1-PBS in a different storage buffer formulation.

E.coli expressed SARS-CoV-2 Nucleocapsid Phosphoprotein (Cat.NO. Ag30676) was subjected to SDS-PAGE followed by western blot with 67666-1-Ig and 67666-2-Ig at various work concentration. This data was developed using the same antibody clone with 67666-1-PBS in a different storage buffer formulation.



Sandwich ELISA standard curve of MP50061-1, Virus 2019-nCoV nucleocapsid phosphoprotein Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67666-1-PBS. Detection antibody: HRP-conjugated 67666-2-PBS. Standard: Ag30676. Range: 0.5-20 ng/mL.