

For Research Use Only

GIMAP2 Monoclonal antibody, PBS Only



Catalog Number: 67720-1-PBS

Basic Information

Catalog Number: 67720-1-PBS	GenBank Accession Number: BC032345	Purification Method: Protein G purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 26157	CloneNo.: 1F2A2
Source: Mouse	UNIPROT ID: Q9UG22	
Isotype: IgG1	Full Name: GTPase, IMAP family member 2	
Immunogen Catalog Number: AG4368	Calculated MW: 337 aa, 38 kDa	
	Observed MW: 33-38 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MOLT-4 cells were subjected to SDS PAGE followed by western blot with 67720-1-Ig (GIMAP2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67720-1-PBS in a different storage buffer formulation.