

For Research Use Only

SAR1B Monoclonal antibody, PBS Only



Catalog Number: 67765-1-PBS

Basic Information

Catalog Number: 67765-1-PBS	GenBank Accession Number: BC002847	Purification Method: Protein G purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 51128	CloneNo.: 1B11A1
Source: Mouse	UNIPROT ID: Q9Y6B6	
Isotype: IgG1	Full Name: SAR1 homolog B (S. cerevisiae)	
Immunogen Catalog Number: AG7718	Calculated MW: 22 kDa	
	Observed MW: 22-25 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, mouse, rat, pig, rabbit

Background Information

SAR1B, also named as SARA2, SARB and GTBPB, belongs to the small GTPase superfamily and SAR1 family. It is involved in transport from the endoplasmic reticulum to the Golgi apparatus. SAR1B is activated by the guanine nucleotide exchange factor PREB. It is involved in the selection of the protein cargo and the assembly of the COPII coat complex.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67765-1-Ig (SAR1B antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67765-1-PBS in a different storage buffer formulation.