

For Research Use Only

# HGS Monoclonal antibody, PBS Only (Capture)



Catalog Number: 67818-4-PBS

## Basic Information

<b>Catalog Number:</b> 67818-4-PBS	<b>GenBank Accession Number:</b> BC003565	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 9146	<b>CloneNo.:</b> 2A9B4
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> O14964	
<b>Isotype:</b> IgG1	<b>Full Name:</b> hepatocyte growth factor-regulated tyrosine kinase substrate	
<b>Immunogen Catalog Number:</b> AG28610	<b>Calculated MW:</b> 86 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Product Information

67818-4-PBS targets HGS as part of a matched antibody pair:

MP50245-3: 67818-4-PBS capture and 67818-3-PBS detection (validated in Cytometric bead array)

MP50245-4: 67818-4-PBS capture and 67818-5-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50245-3, HGS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67818-4-PBS. Detection antibody: 67818-3-PBS. Standard:Ag28610. Range: 3.125-100 ng/mL



Cytometric bead array standard curve of MP50245-4, HGS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67818-4-PBS. Detection antibody: 67818-5-PBS. Standard:Ag28610. Range: 1.563-100 ng/mL