

For Research Use Only

HGS Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67818-5-PBS

Basic Information

Catalog Number: 67818-5-PBS	GenBank Accession Number: BC003565	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 9146	CloneNo.: 1D5E6
Source: Mouse	UNIPROT ID: O14964	
Isotype: IgG2a	Full Name: hepatocyte growth factor-regulated tyrosine kinase substrate	
Immunogen Catalog Number: AG28610	Calculated MW: 86 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

67818-5-PBS targets HGS as part of a matched antibody pair:

MP50245-2: 67818-2-PBS capture and 67818-5-PBS detection (validated in Cytometric bead array)

MP50245-4: 67818-4-PBS capture and 67818-5-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50245-2, HGS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67818-2-PBS. Detection antibody: 67818-5-PBS. Standard:Ag28610. Range: 1.563-100 ng/mL



Cytometric bead array standard curve of MP50245-4, HGS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67818-4-PBS. Detection antibody: 67818-5-PBS. Standard:Ag28610. Range: 1.563-100 ng/mL