

For Research Use Only

PPARGC1B Monoclonal antibody, PBS Only



Catalog Number: 67821-1-PBS

Basic Information

Catalog Number:

67821-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28861

GenBank Accession Number:

BC132971

GeneID (NCBI):

133522

UNIPROT ID:

Q86YN6

Full Name:

peroxisome proliferator-activated receptor gamma, coactivator 1 beta

Calculated MW:

1023 aa, 113 kDa

Observed MW:

150 kDa

Purification Method:

Protein G purification

CloneNo.:

2D2E8

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, pig, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67821-1-Ig (PPARGC1B antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67821-1-PBS in a different storage buffer formulation.



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with 67821-1-Ig (PPARGC1B antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67821-1-PBS in a different storage buffer formulation.