

For Research Use Only

EMCN Monoclonal antibody, PBS Only



Catalog Number: 67854-1-PBS

Basic Information

| | | |
|--|--|---|
| Catalog Number: 67854-1-PBS | GenBank Accession Number: BC017781 | Purification Method: Protein G purification |
| Size: 100ug , Concentration: 1 mg/ml by Nanodrop; | GeneID (NCBI): 51705 | CloneNo.: 2D8G3 |
| Source: Mouse | UNIPROT ID: Q9ULCO | |
| Isotype: IgG1 | Full Name: endomucin | |
| Immunogen Catalog Number: AG31146 | Calculated MW: 19 kDa, 27 kDa | |
| | Observed MW: 120 kDa | |

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, Mouse

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HUVEC cells were subjected to SDS PAGE followed by western blot with 67854-1-Ig (EMCN antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67854-1-PBS in a different storage buffer formulation.