

For Research Use Only

HABP2 Monoclonal antibody, PBS Only



Catalog Number: 67867-1-PBS

Basic Information

Catalog Number: 67867-1-PBS	GenBank Accession Number: BC031412	Purification Method: Protein A purification
Size: 100ug, Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 3026	CloneNo.: 1A7C2
Source: Mouse	UNIPROT ID: Q14520	
Isotype: IgG2b	Full Name: hyaluronan binding protein 2	
Immunogen Catalog Number: AG3981	Calculated MW: 560 aa, 63 kDa	
	Observed MW: 60,63 kDa	

Applications

Tested Applications:
WB, IF, Indirect ELISA

Species Specificity:
Human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human plasma were subjected to SDS PAGE followed by western blot with 67867-1-Ig (HABP2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67867-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using HABP2 antibody (67867-1-Ig, Clone: 1A7C2) at dilution of 1:800 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 67867-1-PBS in a different storage buffer formulation.