

For Research Use Only

WIF1 Monoclonal antibody, PBS Only

Catalog Number: 67904-1-PBS



Basic Information

Catalog Number: 67904-1-PBS	GenBank Accession Number: BC018037	Purification Method: Protein G purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 11197	CloneNo.: 2G6A3
Source: Mouse	UNIPROT ID: Q9Y5W5	
Isotype: IgG1	Full Name: WNT inhibitory factor 1	
Immunogen Catalog Number: AG23972	Calculated MW: 42 kDa	
	Observed MW: 41 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67904-1-Ig (WIF1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Alpha Tubulin Monoclonal antibody (HRP-66031) as loading control. This data was developed using the same antibody clone with 67904-1-PBS in a different storage buffer formulation.