

For Research Use Only

NAT1 Monoclonal antibody, PBS Only



Catalog Number: 67942-1-PBS

Basic Information

Catalog Number: 67942-1-PBS	GenBank Accession Number: BC047666	Purification Method: Protein G purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 9	CloneNo.: 1E1E4
Source: Mouse	UNIPROT ID: P18440	
Isotype: IgG1	Full Name: N-acetyltransferase 1 (arylamine N-acetyltransferase)	
Immunogen Catalog Number: AG6280	Calculated MW: 34 kDa	
	Observed MW: 32 kDa	

Applications

Tested Applications:
WB, Indirect ELISA
Species Specificity:
Human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67942-1-Ig (NAT1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67942-1-PBS in a different storage buffer formulation.